

PHYSICS

CERTIFICATE OF ACHIEVEMENT

THIS CERTIFICATE IS AWARDED TO


FOR SUCCESSFULLY COMPLETION OF
HIS/HER PROJECT WORK.

Affix Your Passport
Size Photo Here

Name of the Student - Rahul
Roll no. 04820902819 Class - ECE (Ist year)

Examination Centre _____


Date of the Practical Examination _____



Signature

Signature

S.No	Experiment Description	Page No.	Experiment Date	Submission Date	Remarks
1	To plot a graph b/w distance of knife edge from centre of gravity.	1	19/08/19		
2	To find wavelength of He-Ne laser source	4	21/08/19		
3	Find dispersive power of glass prism using spectrometer.	6	28/08/19		16/10/19
4	To find n.A of optical fibre.	9	11/09/19		
5	To find wavelength of sodium light by Newton's ring.	11	04/10/19		
6	To det. specific rot. of sugar/mol with half shade polarimeter	14	9/10/19		30/10/2019.
7	To determine planck's constant.	16	16/10/19		
8	To det. wavelength of spectral lines using diff. grating.	18	23/10/19		



Bar pendulum

Aim \Rightarrow To plot a graph b/w the distance of the knife edge from the centre of the gravity and time period of bar pendulum from graph, find.

1. The acceleration due to gravity.
2. The radius of gyration and the moment of inertia of the bar about an axis.

Apparatus :-

A bar pendulum, meter scale, stopwatch, graph paper and sharp edge.

Theory and formula :-

The bar pendulum consist of a uniform metal bar about a meter long with holes drilled along its length at equal interval of about 5cm.

If, T is the time period

L is the length of the pendulum

g is the acceleration due to gravity.

Then, we know that, $T = 2\pi \sqrt{\frac{l}{g}}$

$$T^2 = 4\pi^2 \frac{l}{g}$$

$$g = \frac{4\pi^2 l}{T^2}$$

Radius of gyration, $k = \sqrt{0.125}$

Observations:

Leapt count of the stop watch = 0.1 sec

Leapt count of metre scale = 0.01 m

Mass of the bob pendulum = 1.415 kg

S.NO.	Distance from from O.C (cm)	Time taken for 20 oscillations (sec)	Time period (T) (s)
1	End - A 245 cm	32.5 sec	1.625
2	240 cm	31 sec	1.55
3	235 cm	30.5 sec	1.525
4	230 cm	30.5 sec	1.525
5	225 cm	29 sec	1.45
6	20 cm	3.32 sec	1.56
7	15 cm	3.37 sec	1.687
8	10 cm	3.85 sec	1.925
9	5 cm	5.25 sec	2.625
END - B			
1	5 cm	3.2 sec	2.725
2	10 cm	3.5 sec	1.925
3	15 cm	3.75 sec	1.775
4	20 cm	4.75 sec	2.375
5	25 cm	5 sec	1.6
6	30 cm	5.75 sec	1.475
7	35 cm	6.25 sec	1.5
8	40 cm	7 sec	1.537
9	45 cm	7.5 sec	1.575


where l' and l'' are distance of the points of suspension and point of oscillation.

$$I = mk^2$$

where, m = Mass of bar pendulum
and, k = Radius of gyration.

Procedure:-

1. Put the knife-edge in 1st hole on both sides of the bar.
2. Suspend the pendulum from the knife edge on the side end A of bar.
3. Set the pendulum into oscillations with small amplitude ($\sim 3^\circ - 5^\circ$).
4. Note the time taken for 20 oscillations and measure the distance of the hole from the C.G of the bar.
5. Repeat the observation with knife edges in the 2nd, 3rd and 4th holes.
6. Now suspend the bar from the knife edge on side B.
7. Note the time taken for 20 oscillations with knife edges in the 1st, 2nd, 3rd hole and the distance of edges from the C.G of the bar.
8. Note down the mass of bar pendulum.



Distance from C.G (m)

$$AC = 34 + 23 = 57 \text{ cm} \Rightarrow l' = \frac{57}{2} = 28.5$$

$$BD = 22 + 33 = 55 \text{ cm} \Rightarrow l'' = \frac{55}{2} = 27.5$$

$$l_h = \frac{AC + BD}{2} = \frac{57 + 55}{2} = 56$$

$$g = \frac{4\pi^2 l}{T^2} = \frac{4 \times (3.14)^2 \times 56}{(1.5)^2} = 981.5779 \text{ cm/s}^2 = 9.815779 \text{ m/s}^2$$

$$\% \text{ error} = \left[\frac{9.815779 - 9.81}{9.81} \right] \times 100 = 0.0571\%$$

$$K = \sqrt{l' l''} = \sqrt{0.28 \times 0.27} = 0.274 \text{ m}$$

$$KI = mk^2 = 1.45 \times (0.27)^2 = 0.1057 \text{ kg/m}^2$$


Result

- 1 The acceleration due to gravity, $g = 9.8157 \text{ m/s}^2$
- 2 The radius of gyration = 0.274 m
- 3 Moment of Inertia = 0.1057 kg/m^2

Precautions:-

- 1 The amplitude of vibration should be small.
- 2 The bar pendulum should be hang in vertical plane and knife edge should be in horizontal plane.
- 3 The knife should be sharp.
- 4 The time of oscillations should be counted for atleast 20 oscillations.

8
16/10/2019



Aim \Rightarrow To find the wavelength of He-Ne laser source using transmission diffraction grating.

Apparatus:-

He-Ne laser source, diffraction grating, screen, stands, meter scale and power supply.

Theory and formula:-

When a parallel beam of light of wave length λ be the incident normally on a diffraction grating. The principal maxima of n^{th} order is given by:-

$$(e+d) \sin \theta_n = n\lambda$$

$$\Rightarrow \lambda = \frac{(e+d) \sin \theta_n}{n}$$

where, e = width of each slit,

d = width of each opaque space b/w two slits

θ_n = an angle of diffraction

n = order

$(e+d)$ = grating element (inch)

$$(e+d) = \frac{2.54}{N} \text{ per cm.}$$

Normally, N = No. of ruled lines on a grating
[12,000 to 36,000 Per inch].

Observations

No. of lines on diffraction grating, $n = 1000$ per mm
 No. of lines on diffraction grating, $n = 1000$ per mm
 grating element, $(e+td) = \frac{1000}{1000} = 1$ per mm
 $\therefore e = 1 - td = 1 - 0.001 = 0.999$
 $\therefore e = 0.999$ per mm

Order (n)	Diffraction b/w screen (cm) (D)	Distance of max from central point O (cm) (d)	Distance of max from Sun On (dm)
1	98.28	6.5	0.065989
2	98.07	13.5	0.13637
3	98.08	19.5	0.19500
4	80.99	15.7	0.18030
5	79.57	33.4	0.38712
6	79.67	33.3	0.38911

We know,

$$1 \text{ A}^\circ = 1 \times 10^{-8} \text{ cm}$$

$$\lambda_1 = 6.5989 \times 10^{-5} \text{ cm}$$

$$\lambda_2 = 6.8188 \times 10^{-5} \text{ cm}$$

$$\lambda_3 = 6.83000 \times 10^{-5} \text{ cm}$$

$$\lambda_4 = 6.90300 \times 10^{-5} \text{ cm}$$

$$\lambda_5 = 6.93500 \times 10^{-5} \text{ cm}$$

$$\lambda_6 = 6.29703 \times 10^{-5} \text{ cm}$$

$$\lambda = \lambda_1 + \lambda_2 + \lambda_3 + \lambda_4 + \lambda_5 + \lambda_6$$

$$= 4.17840 \times 10^{-5} \text{ cm}$$

$$= 4178.4 \text{ A}^\circ$$

$$\% \text{ error} (\delta) = (4178.4 - 6328) / 6328$$

$$= -34.0$$

$$\Rightarrow \frac{2152.6}{6328} \times 100 = 34.0$$

Procedure :-

- 1 Set the He-Ne laser source horizontally.
- 2 Place the diffraction grating between He-Ne laser source and screen.
- 3 Find out the principal maximum, 1st and 2nd order maximum.
- 4 Now measure the distance b/w diffraction grating and screen that is D (cm).
- 5 Measure the distance from principal maximum to 1st maximum that is d (n=1).
- 6 Now find out the value of $\sin\theta$, by formula.
- 7 Repeat the steps 1-6 three times for different distance D and note down the reading.


Result :-

- 1 The wave length of He-Ne laser = 6532.5 Å°

Precautions:-

- 1 Never look directly into the laser source.
- 2 The He-Ne laser should be horizontal.
- 3 Adjust the distance b/w diffraction grating and screen properly.
- 4 The diffraction grating should be handled by the edges.

 16/10/2019



Spectrometer

Aim \Rightarrow To determine the dispersive power of a glass prism using spectrometer and mercury source.

Apparatus: \Rightarrow Spectrometer, mercury source, glass prism, magnifying lens and spirit level.

Theory and Formula

It is a measure of the changes of the refractive index of the medium with wave length. The dispersive power of a prism is given by

$$\omega = \frac{m_v - m_x}{m - 1}$$

where, m_v = refractive index for violet colour.

m_x = refractive index for red colour.

m = Mean of refractive index of violet and red colour.

$$\text{i.e. } m = \frac{m_v + m_x}{2}$$

Calculate value of m_x

$$m_x = \frac{\sin A + \sin m}{2}$$

$$\sin A/2$$

where S_m = angle of minimum deviation.

Observation table

S. No.	Vernier Reading	Reading focused on left faced (θ_1)	Reading focused on right faced (θ_2)	$\theta = \theta_1 - \theta_2$	$A = 0/$
1	v_1	$33^{\circ} 6'$	$272^{\circ} 5'$	$121^{\circ} 1'$	$60^{\circ} 31'$
	v_2	$213^{\circ} 21'$	$92^{\circ} 41'$	$120^{\circ} 58'$	$60^{\circ} 29'$

Mean $A = 60^{\circ} 30'$ (approx)
 $= 60.5^{\circ}$

1) For Violet light

S. No.	Vernier Reading	Telescope reading (min. dev. pos. A)	Direct (B) reading	Difference $D_m = A - B$
1	v_1	$210^{\circ} 9'$	$349^{\circ} 6'$	$51^{\circ} 3'$
	v_2	$220^{\circ} 3'$	$169^{\circ} 5'$	$50^{\circ} 58'$

Mean $D_m = 51^{\circ} 34'$

2) For Red light

S. No.	Vernier Reading	Telescope reading min. dev. pos. A	Direct reading (B)	Difference $D_m = A - B$
1	v_1	$402^{\circ} 8'$	357.9	$44.9'$
	v_2	$223^{\circ} 1'$	178.7	$47.4'$

Mean $D_m = 46.15'$

Procedure:-

1 For calculation of A.

- (1) Adjust the spectrometer to get parallel light beam and check prism table is horizontal to spirit level.
- (2) Now put prism on prism table such that two refracting edges should point towards the light source.
- (3) Rotate the telescope in the direction of and then along 90° to see the diffraction pattern that is the slit of white light on the crosswire of the telescope. Write down the position of these two images from the vernier scales. The difference b/w these two points gives twice the value of A.

2 For measurement of δ_m

- (1) Now, remove the prism from the table. Bring the telescope to the straight position and note down the reading.
- (2) Choose any two colours in the spectrum. Rotate the telescope to right side and look for the slit of the two colour of the crosswire. Note down their positions from the vernier scale.
- (3) The d/f b/w position of a particular colour with position of straight view will give the angle of min. deviation for the particular colour.
- (4) Now calculate u and w .

$$\mu_v = \frac{\sin(60^\circ 5' + 51^\circ 34')}{2} = \frac{\sin 86^\circ 02'}{2}$$

$$\sin(60^\circ 5'/2) = \frac{\sin 30^\circ 25'}{2}$$

$$= \frac{0.811}{2} = 0.507$$

$$\therefore \mu_v = 1.600$$

$$\mu_x = \frac{\sin(60^\circ 5' + 46.15)}{2}$$

$$\sin 53^\circ 32'$$

$$\sin(60^\circ 5'/2) = \frac{\sin 30^\circ 25'}{2}$$

$$= \frac{0.811}{2} = 0.507$$

$$= 0.507$$

$$= 1.577$$

$$\therefore \mu = \frac{\mu_v + \mu_x}{2} = \frac{1.6 + 1.577}{2} = 1.5885$$

$$= \frac{3.17}{2} = 1.5885$$

$$\omega = \frac{\mu_v - \mu_x}{\mu - 1} = \frac{1.6 - 1.577}{1.5885 - 1} = \frac{0.03}{0.59}$$

$$= 0.05$$

Result:-

The dispersive power of prism = 0.05

Precautions:-

1 The reading should be noted properly.


2 Cross wire should be focussed properly.

3 Prism table should be levelled.

4 Pattern should be narrow and sharp.


Diagram

(1)




Numerical aperture = Apperture angle \times $\sin \theta$

(2)



(3)



N.A. of fibre

Aim \Rightarrow To measure the numerical aperture (NA) of an optical fibre.

Apparatus \Rightarrow Optical fibre, fibre optic chuck, He-Ne laser source, microscopic objective, screen, graph paper, scale.

Theory and formula

Snell's law can be used to calculate the maximum angle with which light will be accepted into and conducted through a fiber

$$N.A = \sin \theta_n = (n_1^2 - n_2^2)^{1/2}$$

where, $\sin \theta_n$ is numerical aperture
 n_1 and n_2 are the refractive indices of the core and the cladding.

Procedure:-

- 1 Mount both the ends of the optical fibre on the fibre optic chucks
- 2 Couple the light from the He-Ne laser source onto one of the fiber ends using a 10x microscopic objective.

Observations:-

S.NO.	Distance (L) in cm.	Diameter (d) in cm
1	0.7	1.6 cm
2	1.3	2.0 cm
3	1.7	2.4 cm
4	2.3	3.2 cm
5	2.8	3.6 cm
6	3.1	4.0 cm

Calculation:-

$$\text{Mean of distance (L)} = \frac{0.7 + 1.3 + 1.7 + 2.3 + 2.8 + 3.1}{6} = 1.98 \text{ cm}$$

$$\Rightarrow \frac{11.8}{6} = 1.98 \text{ cm}$$

$$\text{Mean of diameter (d)} = \frac{1.6 + 2.0 + 2.4 + 3.2 + 3.6 + 4.0}{6} = \frac{16.8}{6} = 2.80 \text{ cm}$$

$$(180^\circ / N \cdot A) = \sin^{-1} \left[\frac{2}{2 \times L} \right]$$

$$= \sin^{-1} \left[\frac{2.8}{2 \times 1.98} \right] = \sin^{-1} \left[\frac{1.4}{3.96} \right]$$

$$= \sin^{-1} (35.26) = 0.577$$

- 3 Place the screen at some distance from the output end of the fiber such that it is \perp to the axis of fiber.
- 4 Now move screen towards output end of a optical fibre, such that a circular spot is formed on the screen.
- 5 Measure the distance b/w the output end of optical fibre and screen. Let this be L . Also measure diameter of circular spot formed on screen.
- 6 Use the formula, $NA = \sin \alpha [\tan^{-1} (\lambda / 2L)]$
- 7 Repeat the above procedure for different values of L and λ and calculate the average of NA .


Result:-

N.A of given optical fibre is 0.577

Precautions:-

- 1 Mounting and coupling should be carefully done.
- 2 The laser source shouldn't fall into eye directly.
- 3 The screen should be held vertically.
- 4 As far as possible, experiment should be conducted in dark room environment.

16/10/19



Aim \rightarrow To determine the wavelength of sodium light by Newton's ring.

Apparatus :-

A plane convex lens of large radius of curvature, a travelling microscope, a sodium lamp, plane glass plate, spectrometer and reading lens.

Theory and formula :-

Let R be the radius of curvature of the surface in contact with the plate, λ be the wavelength of light used. D_n and D_{m+n} the diameter of n^{th} and $(m+n)^{\text{th}}$ dark rings respectively, then

$$D_n^2 = 4n\lambda R$$

and $D_{m+n}^2 = 4(m+n)\lambda R$

$$\text{or } D_{m+n}^2 - D_n^2 = 4m\lambda R$$

$$\Rightarrow \lambda = \frac{D_{m+n}^2 - D_n^2}{4mR} \quad \textcircled{1}$$

Using eq. ①, λ can be determined.

Procedure :-

- 1 All the eyepiece of microscope adjusted on its cross wires.
- 2 Now, distance of microscope from the film is adjusted at the rings with dark centre in well focus.

Observations:

S. NO.	Microscope reading LHS (cm)	Microscope reading RHS (cm)	D (cm)	$D^2(a+b)^2$	$D_{n+p}^2 - D_n^2$	$P = 1 \text{ cm}$ and Mean
1	45.705	40.4	5.303	28.143	6.669	60.099
2	45.1	40.703	4.146	16.948	20.252	44.4215
3	46.103	40.509	5.594	31.292	9.616	44.4215
4	46.801	41.109	5.612	32.464	21.149	44.4215
5	46.908	41.007	5.400	31.812	3.812	44.4215
6	46.906	40.807	6.009	37.107	10.107	44.4215
7	47.303	40.907	6.396	40.908	10.908	44.4215
8	47.705	41.106	6.599	43.553	13.553	44.4215

Calculation:

$$\lambda = \frac{(D_{n+p}^2 - D_n^2)}{4PR} \text{ mean} = \frac{14.4215}{4 \times 4 \times 2} \times 10^{-3}$$

$$= 0.4506718 \times 10^{-3} \text{ mm} = 4506.718 \text{ Å}^2 \quad [1 \times 10^{-7} \text{ mm} = 1 \text{ Å}]$$

But wavelength of sodium light = 589.3 nm
 Or 5890 Å

$$\% \text{ error} = \left[\frac{5890 - 4506.718}{5890} \right] \times 100 = 23.49 \%$$

- 3 The centre of cross wires is adjusted at the centre of fringe pattern.
- 4 By counting no. of fringes, the microscope is moved to the extreme left of the pattern and the cross wire is adjusted tangentially in the middle of n^{th} (say 20th) bright or dark fringes.
- 5 The reading of micrometer screw is noted.
- 6 The microscope is now moved to the right and readings of micrometer screw are noted successively at $(n-2)^{\text{th}}$ (say 18th), $(n-4)^{\text{th}}$ (say 16th), etc.
- 7 Now, plotted a graph b/w a no. of rings (n) and square of corresponding diameter.

From the graph, $\frac{D_{m+n}^2 - D_n^2}{m} = \frac{AB}{CD} \quad \text{--- } ②$


Radius of planoconvex lens, (R) = $\frac{l^2}{6n} + \frac{b}{2} \quad \text{--- } ③$

Result

The wavelength of sodium light by Newton's rings experiment is 3725 Å .

Precautions:-

- 1 Glass plate and lens should be cleaned thoroughly.




16/10/19

- 2 The lens used should be large radius of curv.
- 3 Source of light used should be extended one.
- 4 Before measuring the diameter of the rings, the range of the microscope should be properly adjusted.
- 5 Cross wire should be focused on a bright ring tangentially.
- 6 Radius of curvature should be measured accurately.

Ansatz der optischen Achse nach dem
Sagittalschnitt bei einem Kreiselliptisch

Diagramm: Einstrahlung parallel zum Spalt des
Linsenobjektivs und Linsenachse senkrecht zu
dieser. Durch ein Polarisator wird die
Polarisation des Lichtes so geändert, dass
die Linsenachse und die Blendenachse (die Längs-
achse) übereinstimmen.



half shade
Polarimeter

Aim \Rightarrow To determine the specific rotation of cane sugar salt with the help of a half shade polarimeter.

Apparatus:-

Polarimeter, sodium light, flask, beaker, funnel, pipette, sugar and distilled water.

Formula/Theory

$$\text{The specific rotation of sugar salt, } S = \frac{\Theta V}{l m}$$

where, Θ - rotation produces in degrees,

m = mass of the sugar dissolved,

l = length of tube in decimeter.

V = Volume of the soln.

Procedure:-

- 1 Fill the polarimeter tube with water minimizing the air gap and place it in proper posn.
- 2 The slit is illuminated by sodium light and focus the eye piece such that the field of view is sharp and unequal illuminated columns are seen in the two halves of the field of view. The analyser is rotated till the two position or colours change into homogeneous colour.
- 3

Observation:

Analyser reading with pure water.		Mean(y) $\frac{x+y}{2}$	Angle of ref. on gm/cc	Analyser reading with sugar soln.	
Clockwise	Anticlock	deg.	deg.	Clockwise	Anticlock
Total(x) = $m \cdot s + v \cdot s \times LC$	Total(y) $= m \cdot s + v \cdot s \times LC$			Total(x') = $m \cdot s + v \cdot s \times LC$	Total(y') $= m \cdot s + v \cdot s \times LC$
1 103.7	201.2	152.45	$\frac{5}{100} = 0.05$	138.1	(44.1)
2 198.4	223.3	211.35	$\frac{2.5}{100} = 0.025$	125.5	170.9
3 138.9	232.4	185.15	$\frac{1.25}{100} = 0.0125$	100.40	286.2

Calculation

Least count of polarimeter = 0.1 deg.

Mass of sugar (m) = 10 g

Volume of water (v) = 100 ml

length of polarimeter tube (l) = 2 dm

$$S = \frac{0.1 \times 2}{100} = \frac{0.1 \times 2}{100} = 0.002 \text{ dm/gm/cm}^3$$

$$\% \text{ error} = \left[\frac{\text{Standard value} - \text{exp. value}}{\text{Standard value}} \right] \times 100$$

$$= \left[\frac{66 - 56.75}{66} \right] \times 100 = 14.01\%$$

- 4 Take the measurement on circular and Vernier scales. Rotate the analyser about 180° for some colour and take the measurement again.
- 5 Prepare the sugar soln of known concentration (say 10 gm in 100 ml). Fill the polarimeter tube with sugar solution.
- 6 The difference in the observations taken in distilled water and in sugar soln is the angle of rotation (θ) for one concentration.
- 7 Repeat the experiment by changing concentration of sugar soln and take atleast five sets of readings.
- 8 Plot the graph b/w concentration and the angle of rotation. It will be a straight line.


Result:-

The specific rotation of cane sugar is 56.25° /dm³ cm/cm²

Precautions:-

- 1 There should not be any air bubbles in polarimeter.
- 2 Caps must be put on the tube tightly without undue strain.
- 3 Water used should pure/ distilled.
- 4 Sugar must be dust free.
- 5 Before filling the tube, it should be rinsed with the soln to be filled.

Diagram:-



Observations:-

S. No.	Filter Colour	Wavelength (λ) nm	Frequency ($V = f$ (THz))	Stopping Potential (E)
1	Blue	460	625.47	1.16
2	Orange	500	600	0.89
3	Green	540	555.5	0.73
4	Yellow	570	526.3	0.55
5	Red	635	472.4	0.34

Calculation:-

According to Einstein's eqn:-

$$h\nu = h\nu_0 + e\Delta E_0$$

$$h(\nu - \nu_0) = e\Delta E_0 \Rightarrow h = \frac{e\Delta E_0}{\nu - \nu_0} = e \tan \theta$$

Aim \Rightarrow To determine Planck's constant.

Apparatus:-

A vacuum photo cell, mercury lamp, five optical fibres (Red, yellow, green, blue, orange).

Theory

The experiment is based on the phenomenon of photo-electric effect. A/c to this when a radiation of energy $h\nu$ incident on metal surface then a photoelectron is emitted. A/c to Einstein eqⁿ.

$$h\nu = \frac{1}{2}mv_m^2 + W_0$$

$$eE_0 = \frac{1}{2}mv_m^2 \\ = h\nu - h\nu_0$$


$$\Rightarrow E_0 = \frac{h\nu}{e} - \frac{h\nu_0}{e}$$

If θ is inclination of line from frequency axis, then,

$$\tan\theta = \frac{h}{e} = \frac{\Delta E_0}{\Delta V}$$

$$\Rightarrow h = \frac{e\Delta E_0}{\Delta V}$$

Graph:



$$\therefore \text{Slope} = \frac{y_2 - y_1}{x_2 - x_1}$$

$$= \frac{0.55 - 0.34}{(526.3 - 472.9) \times 10^{-2}} = \frac{0.21}{53.9 \times 10^{-2}}$$

$$= 3.89 \times 10^{-15} \text{ V/Hz}$$

$$S_e, S_g, h = \text{Slope} \times e$$

$$= 3.89 \times 10^{-15} \times 1.6 \times 10^{-19}$$

$$= 6.224 \times 10^{-34} \text{ JS}$$

$$\% \text{ error} = \left[\frac{6.626 - 6.224}{6.626} \right] \times 100 = 6.06\%$$

Procedure:-

- 1 Note the wavelength of given filter.
- 2 Switch on the Planck's constant measuring kit
- 3 Keep a fix distance b/w source and cathode.
- 4 Place the filter towards cathode and vary the anode voltage negatively starting from zero. Measure the corresponding voltage and current till the current is equal to zero. The voltage at zero current is called as stopping potential.
- 5 Repeat step 4, using different filters.
- 6 Make a table between frequency of light emitted by filter and stopping potential.
- 7 Make a graph: b/w E and V and find the slope.

Result:-


The value of planck's const. using photocell is 6.224×10^{-34} J S

Standard value of Planck's const. = 6.626×10^{-34} JS

Precautions:-

- 1 Avoid the extra light falling on the photocell.
- 2 Set a experiment at a place where it can be avoided from outside disturbances.
- 3 Phototube shouldn't be exposed to direct sunlight.
- 4 Power should switch off on finishing the experiment.

✓
30/10/2019.



Order of spectrum	Colour of spectrum	Kind of Venuer	LHS Reading of spectrometer			RHS Reading of spectrometer			$\frac{a-b}{2}$ degree
			MS Reading	VS Reading	Total degree	MS Reading	VS Reading	Total degree	
1st	Violat	V ₁	280	26	200.43	306.5	32	307.03	26.60
		V ₂	100	17	106.28	176.5	33	124.05	26.77
Order	Green	V ₁	276.3	92	276.51	310	26	310.43	33.92
		V ₂	96.5	14	97.18	130	23	130.38	32.58
Yellow		V ₁	275.5	13	275.71	311	29	311.48	35.77
		V ₂	95.5	9	95.65	131	30	131.50	35.85

Aim \rightarrow To determine the wavelength of spectral lines using diffraction grating.

Apparatus:-

Spectrometer, mercury lamp, prism, diffraction grating, reading lens.

Formula / Theory

The wavelength of any spectral line,

$$\lambda = \frac{(a+b) \sin \theta}{n}$$

where, $(a+b)$ = grating element

θ = angle of diffraction

n = order of spectrum

Procedure:-

- 1 Prism table should be leveled.
- 2 Telescope and collimator are adjusted for parallel rays.
- 3 The axis of the telescope and collimator must intersect the vertical axis of rotation of the prism table.
- 4 Grating should be normal to the axis of collimator.

Calculation:-

No. of rulings on the grating = $N = 15000$

Least count of spectrometer = $\frac{1}{N} = \frac{1}{15000} \text{ of } 60^\circ = 0.00067^\circ$

Grating element ($a+b$) = $\frac{2.54}{N} = \frac{2.54}{15000} \text{ cm}$

$= 1.693 \times 10^{-6} \text{ cm}$

where, $N = \text{no. of lines per inch}$

now, for 1st order $\lambda_1 = (a+b) \sin \theta$

For violet, $\lambda_V = 1.693 \times 10^{-6} \sin(26.68)$

$$= 1.693 \times 10^{-6} \times 0.449$$

$$= 0.7601 \times 10^{-6} \text{ cm}$$

$$= 0.7601 \times 10^{-2} \text{ Å}$$

For green, $\lambda_G = 1.693 \times 10^{-6} \sin(33.25)$

$$= 1.693 \times 10^{-6} \times 0.548$$

$$= 0.927 \times 10^{-6} \text{ cm}$$

For yellow,

$$\lambda_Y = 1.693 \times 10^{-6} \sin(35.81)$$

$$= 1.693 \times 10^{-6} \times 0.5850$$

$$= 0.990 \times 10^{-6} \text{ cm}$$

$$= 0.990 \times 10^{-2} \text{ Å}$$

- 5 After obtaining spectrum, rotate the telescope to the left side of the direct image.
- 6 Now, adjust the cross wire on violet, green, and red spectral lines turn for the 1st order.
- 7 For each setting, the reading of Vernier's should be noted down.
- 8 Now, rotate the telescope to the right side of the direct image.
- 9 Again, note down the position of both the verniers turn by turn for different colours for 1st order.
- 10 Find out the difference of the same kind of verniers (v_1 from v and v_2 from v_2) for each spectral lines.
- 11 Repeat all the observations for second order in similar way.

Result :-

The wavelength of various colours (violet, green, yellow) of white light (mercury light) with the help of plane transmission grating are $0.760 \times 10^{-2} \text{ A}^\circ$, $0.927 \times 10^{-2} \text{ A}^\circ$ and $0.990 \times 10^{-2} \text{ A}^\circ$ respectively.

Precautions:-

- 1 Rack and pinion arrangement of telescope once set for parallel rays should not be disturbed throughout out the experiment.

