

13/01/2021

DELTA Pg No.

Date / /

## EXPERIMENT - I

**Objective** - To plot a graph between distance of the knife edge from the centre of gravity and time period of bar pendulum. From the graph find (a) the acceleration due to gravity  
(b) the radius of gyration.

**Apparatus** - Bar pendulum, a small metal wedge, a spirit-level, a telescope, a stop watch, a metry rod, a rigid support fixed on the wall.

**Theory** - A bar pendulum consists of a uniform metal-bar having holes drilled along its length symmetrically on either side of the centre of gravity at intervals of about 5 cm. Two knife-edges are placed symmetrically with respect to C.G. at A and B. The bar is allowed to oscillate about a horizontal knife-edge passing through each of holes in turn beginning with the nearest to one-end and the time period is noted in each case. If a graph is plotted between the distance of the knife edge from the centre of gravity and the corresponding time period T, it has the shape shown in the figure.

It has a horizontal line ABCDE is drawn, at which the graph is present ABD and E about which the time period is the same for A as the centre of suspension D is the centre of oscillation and for E as the centre of suspension, B is the centre of oscillation.

Stand to hold bar pendulum.




Figure .


The distance AD and BE then give the length of its equivalent simple pendulum. If T is the corresponding time period,  $T$ ,

$$T = 2\pi \sqrt{\frac{L}{g}} \quad T^2 = \frac{4\pi^2 L}{g}$$

The radius of gyration,  $k = \sqrt{LL'}$ ,  
where  $L = AD + \frac{BE}{2}$  and  $L' = BE + \frac{CD}{2}$

### Procedure

1. Shoot/fix/split/level Balance the bar on the sharp wedge mark it C at A.
2. Fix the knife edge in the outermost holes on either end of the bar parallel to each other with their sharp edge pointing towards A.
3. Check the spirit level that the glass plates & fixed on the rigid bracket on plate A fixed on the rigid balance on the wall are horizontal. Suspend the pendulum by placing the knife-edge on the side A perpendicular to the slot of the bracket and to that pendulum, hang vertically.
4. Adjust the eye-piece of the telescope so that the wall-wire are clearly visible through it. Focus the telescope on the bar from about one meter distance and see that the point of intersection of the concave concides with the line AB or A



Calculations from ABCD.

(A) Acceleration due to gravity

$$L = \frac{AO+BE}{2} = \frac{0.60 + 0.60}{2} = 0.60 \text{ m},$$

$$T = 1.63 \text{ sec}$$

$$g = \frac{4\pi^2 L}{T^2} = \frac{4 \times (3.14)^2 \times 0.60}{(1.63)^2} = 9.6 \text{ ms}^{-2}$$

(B) Radius of gyration  $R_g$

$$R_g = \sqrt{\frac{AO+BE}{2}} = \sqrt{0.60 + 0.60} = 0.60 \text{ m}$$

$$J' = \frac{BC+CD}{2} = 0.70 \text{ m}^2$$

$$R = \sqrt{J'/I} = \sqrt{(0.15)(0.70)} = \sqrt{0.08} = 0.28 \text{ m}$$

Result :- 1. The acceleration due to gravity  
 $g = 9.6 \text{ ms}^{-2}$

2. Radius of gyration,  $R = 0.28 \text{ m}$

(A) From  $A = I = Cg^2$   
 Acceleration due to gravity  
 $G = A'D' + B'C' = 0.64 + 0.64 = 0.64$   
 $T = 1.63 \text{ sec}$   
 $g = \frac{4\pi^2 L}{T^2} = \frac{4 \times (3.14)^2 \times 0.60}{(1.63)^2} = 9.5 \text{ ms}^{-2}$

(B) Radius of gyration  $R_g$   
 $J = A'D' + C'C' = 0.40 + 0.40$   
 $J' = \frac{BC + CD}{2} = 0.15 + 0.15$   
 $R = \sqrt{J'/I} = \sqrt{\frac{0.15 \text{ m}^2}{(0.40)(0.15)}} = \frac{0.15 \text{ m}}{\sqrt{0.04}} = 0.27 \text{ m}$

reference line drawn on the pendulum

5. Displace the bob slightly and release it. It will begin to vibrate. Note the time for one vibration and also measure the length from the point upto the bottom of the first hole i.e. upto the point about which the pendulum turns. Distance can also be measured from one end of the bar.
6. Support the pendulum on the knife-edge B side of the pendulum and take the observation.
7. Now invert the knife-edge in hole 1, 2, 3 on either side of the Cu and note the time for 20 vibrations in each case and measure the corresponding distance from a.
8. Plot a curve between distance from C.A and the corresponding time period T. (Time for one vibration)
9. Draw two horizontal line ABCD and A'B'C'D' and not calculation.
10. Find the max of pendulum to determine the Inertia.

Precaution and sources of errors,

1. The knife-edges should be horizontal and the bob pendulum should hang vertically.
2. Amplitude of vibration should be small.
3. The time of oscillation should be noted after 50 oscillations so that any irregularities of motion will be avoided.
4. Time should be noted with the help of a good stop watch.

Observation Table:

	No -	Distance from (in cm)	Time for 20 vibrations / sec.	Time period
for holes on side	1	45 cm	32.29	1.61
A.	2	40 cm	31.2	1.56
	3	35 cm	30.94	1.54
	4	30 cm	30.94	1.52
	5	25 cm	30.57	1.57
	6	20 cm	30.44	1.64
	7	15 cm	32.55	1.89
	8	10 cm	37.97	
	9	5 cm	37.97	2.4
			48	
for holes on side	1	45 cm	32.20	1.61
B	2	40 cm	31.1	1.56
	3	35 cm	30.96	1.54
	4	30 cm	30.60	1.52
	5	25 cm	30.40	1.57
	6	20 cm	31.51	1.52
	7	15 cm	32.82	1.64
	8	10 cm	37.96	1.89
	9	5 cm	48.02	2.4

Beat Count of stop watch: 0.01 sec.


wing a telescope

5. The knife edges should be sharp
6. The support shouldn't be pecky
7. Smooth and good-sized curve should be drawn

weak-points - The mount may be slightly inaccurate  
as a result of:

- (i) The curvature of the knife edges
- (ii) Finite amplitude of the wing
- (iii) yielding of the support and
- (iv) Resistance of air

20/01/21



## EXPERIMENT - 2

Aim - To determine the Planck's constant.

Apparatus - A vacuum type photoelectric-cell mounted inside a wooden box with a wide opening on the side opposite to the cathode, volt DC supply, a voltmeter (0-3 volt range) & ammeter, a sensitive moving coil galvanometer with lamp and scale arrangement, a key, a mercury lamp and few light filters.


Theory - When a photon of energy  $\nu$  is incident on the emissive surface of the cathode, almost all of its energy is transferred to the electron inside the metal. If this energy is greater than threshold energy  $\omega_0$ , the electron ~~move~~ is emitted.  $\omega_0$  is called the work function of the metal. Above the threshold energy frequency, corresponding to  $\omega_0$ , photoelectrons have a range of energies from 0 to a certain maximum value and this maximum increasing frequency. As  $\omega_0$  remains constant for a given photoelectric-cell,  $E_{max}$  varies linearly with frequency. In this experiment to find out the maximum kinetic-energy  $E_{max}$  at the emitted electrons reverse bias the photoelectric-cell so anode is made negative. It therefore repels the emitted electrons and the current decreases. The negative potential on the anode

\* CALCULATIONS : The graph between the frequency  $\nu$ , along X-axis and the stopping potential  $V_s$  along Y-axis is a straight line.

By formula

$$V_s = \frac{h}{e} \nu - \frac{W_0}{e}$$

↓  
Hence  $R = \frac{h}{e}$



Hence,

Planck's constant ( $h$ ) ;

$$h = \text{slope} \times e$$

$$h = 1.6 \times 10^{-19} \times 0.26 \times 10^{-14}$$

$$h = 4.16 \times 10^{-34} \text{ Js}$$

Now,

Calculation of work function

$W_0$  can be

calculated with

solving we did;

Solving them

$$W_0 = e (\text{slope} \times 2) - eV_0 = e [(0.26 \times 10^{-14} \times 6.34 \times 10^4) - 0.185] = 0.943 \text{ eV}$$

### Result

- Planck's constant -  $4.16 \times 10^{-34} \text{ Js}$

- Actual value of Planck's constant -  $3.11 \times 10^{-34} \text{ Js}$

- Work function of given photoelectric cell -  $0.943 \text{ eV}$

- Actual value of work function of given photocell -  $1.6 \text{ eV}$

is usually measured P.D. with stopping potential  $V_s$  is reached, when the current stops.

Formula used -  $h\nu = \frac{1}{2}mv^2 + w_0$ .

$$E_{max} + w_0$$

This is called Einstein's ~~approximate~~ photoelectric equation when stopping potential is reached:-

$eV_s = E_{max}$ .  $V_s$  is called the stopping potential or the cut off potential because it just stops all electrons from leaving the surface.

$$h\nu = eV_s + w_0 \quad \text{or} \quad eV_s = h\nu - w_0$$

$$eV_s = \frac{h}{c} V - \frac{w_0}{c}$$

Thus, if we draw a graph with frequency  $\nu$  along the x-axis, it would be a straight line with slope equal to  $h/e$  and negative y-intercept on y-axis equal to  $w_0/e$ .

### Procedure.

- ① Make the connections as shown in figure. C and A are respectively the cathode and anode of the photocell.
2. Arrange the mercury lamp just in front of the photocell and set the galvanometer so that the spot moves truly.

- ③ With no light falling on the photocell, insert the key in R, adjust the voltage in the voltmeter to zero and adjust spot in galvanometer to zero position. Remove the key.
- ④ Fit the violet filter in front of the photocell.
- ⑤ Increase the potential on anode slowly by the rheostat. The deflection in the galvanometer decreases. Go on increase the negative potential applied to anode till the spot of light in the galvanometer comes back to zero position.
- ⑥ Repeat the process and take three readings with the same filter.
- ⑦ Change the filter one by one from violet to red and repeat the experiment.
- ⑧ Make a graph with frequency of the filter v along X-axis and stopping potential along y-axis.

#### Precaution and sources of error

1. Care should be taken to ensure that when no light is falling on the cathode, the deflection in the galvanometer is zero.
2. The position of the mercury lamp should not be changed during the experiment.
3. Care should be taken to note down the stopping potential. The voltage across the photocell should be increased very slowly and three readings

should be taken. for each filter.


### Weak - Points

A sharp cut-off magnet be obtained due to Penetration of the residual air in the cell and due to secondary emission and photoelectric emission from anode A.

Aim:

To find the wavelength of He-Ne LASER source using transmission diffraction grating.

Figure.



27/10/21

Date \_\_\_\_\_  
Page \_\_\_\_\_

## EXPERIMENT - 3

Aim:

To find the wavelength of He-Ne LASER source using transmission diffraction grating

Apparatus:

He-Ne Laser, Diffraction, optical bench with four stands  
a screen with a graph paper pasted on it, metre scale  
and a travelling microscope

Formula used:

when a monochromatic beam of wavelength  $\lambda$  is diffracted by a diffraction grating, the  $n$ th order principle maxima is formed at an angle  $\theta_n$  given by  
 $(a+b) \sin \theta_n = n\lambda$

where  $(a+b)$  is the grating constant and

$$\lambda = \frac{(a+b)}{n} \sin \theta_n$$

when laser light is incident on the diffraction grating and diffraction pattern is obtained on the screen formed by a graph paper,  $\sin \theta_n$  can be obtained for different order maxima.

From fig. we have

$$\sin \theta_n = \frac{u_n}{\sqrt{L^2 + u_n^2}} = \frac{u_n}{L} \text{ if } L \gg u_n$$

## Observations

Distance of screen from slit gratings "L" = 61.7 cm.

at  $N = 100 \text{ lines/cm}$

Order of Diffraction (n)	Distance of Minima from central spot (cm) (cm)	$\sin \theta_n$ $n = \frac{m}{d^2 - m^2}$	X (cm)
1.	4.8	0.291	7.700
2.	9.8	0.256	7.800
3.	14.5	0.228	7.600
4.	19.6	0.302	7.550

at  $N = 300 \text{ lines/cm}$

Order of Diffraction (n)	Distance of Minima from central spot (cm) (cm)	$\sin \theta_n$ $n = \frac{m}{d^2 - m^2}$	X (cm)
1	14	0.221	7.300
2	29.5	0.431	7.190
3	59	0.537	7.070

$$\text{Mean} = 7458 \text{ A}^\circ$$

Result:

The wavelength of He-Ne Laser light  $\lambda = 7458 \text{ A}^\circ$

Actual Value =  $6328 \text{ A}^\circ$

% Error = 12%

### Precautions and sources of error:

1. Do not see the laser source directly as it may cause damage to eye
2. The graph used at the screen should be vertically straight
3. The diffraction grating should be very near to the laser source


### Result:

The wavelength of He-Ne laser light  $\lambda = 7458 \text{ Å}^\circ$

Actual value =  $6328 \text{ Å}^\circ$

% error =  $17\%$

Aim: To determine the wavelength of sodium light by Newton's rings



Observations table

- Measurement of diameter of rings,  $D_n$  and  $D_n^2$
- Vernier constant reading of microscope 0.01 mm.

S.No.	Number of bright fringes	Microscope reading with cross wire		Diameter of ring $n-y = D_n$	$(D_n)^2$ mm <sup>2</sup>
		On right u (mm)	On right y (mm)		
1	2	42.54	46.20	3.66	13.40
2	4	42.03	46.71	4.68	21.9
3	6	41.78	47.23	5.45	29.70
4	8	41.50	47.60	6.10	37.21
5	10	41.19	47.95	6.76	45.69
6.	12	40.90	48.25	7.35	54.02

EXPERIMENT - 4

Objective - To determine the wavelength of sodium light using Newton's ring.

Apparatus - Newton's ring apparatus (i.e. plano-concave lens) and an optically plane glass plate placed on a wooden stand having a plane glass plate, a travelling microscope, sodium lamp, a convex lens and a photometer.

Theory - Newton's ring are formed as a result of interference between the wave reflected from the top and bottom of air-film surface forward between the lower convex surface and top of the plane P.

By proper analysis we will have following result

$$\lambda = \frac{D_{nem}^2 - D_n^2}{4mR}$$

where,  $D_{nem}$  is the diameter of  $(n_{em})^{th}$  bright ring,  $D_n$  is diameter of  $n^{th}$  bright ring,  $R$  is radius of curvature of the curved surface of plano-concave lens.

- The radius of curvature  $R$  is correlated with the

★ Calculation:


- The radius of the convex surface of the plano-concave lens is 70 cm

- Wavelength by this method calculated equal to  

$$\lambda = \frac{D_m^2 - D_n^2}{4m}$$
 also  $\tan \theta = \frac{(D_m^2 + D_n^2)}{m}$

$$\therefore \lambda > \frac{\tan 22^\circ}{4R} = \frac{0.404 \times 10^{-6} \text{ m}}{4 \times 170 \times 10^{-3} \text{ m}} \approx 0.0005942 \times 10^{-3} = 594.2 \times 10^{-9} \text{ m}$$

$$\boxed{\lambda = 594.2 \text{ nm}}$$



- Result
- Wavelength of medium light obtained  $594.2 \text{ nm}$
  - Actual wavelength of medium light is  $589 \text{ nm}$
  - Per centage error =  $| \frac{594.2 - 589}{589} | \times 100 = 0.88\%$

With help of spherometer, we're using the formula as

$$R = \frac{l^2 + h}{6l}$$

where  $l$  is the distance between the outer leg of the spherometer &  $h$  is the elevation of tip of the central leg as compared to the outer legs when placed on the convex lens.

#### Precaution and source of error:


- ① The surface of plano-convex lens and plane glass plate should be optically spherical and plane respectively.
- ② The plano-convex lens and the glass plate should be cleaned.
- ③ The convex surface of the plano-convex lens should be large to aperture.
- ④ The sodium light should have an extended source.
- ⑤ Before starting the measurement for diameter of the lens, slow motion of the micrometer only in one direction should be ensured for whole range of observation only to avoid back lash error.
- ⑥ Cross-wire should be reflected on the bright defring tangentially.
- ⑦ The value of ' $h$ ' in the measurement of radius will be measured as

**Result:** • Wave length of sodium light obtained : 594.2 nm

$$\text{• Percentage error: } \left| \frac{594.2 - 589}{589} \right| \times 100 = \boxed{0.88\%}$$

Objective To determine the numerical aperture (NA) of an optical fibre.

Setup



Observation

S.NO	Distance between output end of optical fibre and screen L (mm)	Diameter of circular spot D (mm)	$NA = \frac{D}{\sqrt{D^2 + 4L^2}}$
1	1.18	4	0.8613
2	2.66	8	0.8327
3	5.00	12	0.7682
4	6.80	16	0.7619
5	8.70	20	0.7544

## EXPERIMENT - 05

**Objectives** To determine the numerical aperture (NA) of an optical fibre.

**Apparatus** - Optical fibre, two fibre optic checkers, He-Ne laser, a microscopic objective and a screen with a graph paper pasted on it.

**Theory** - Numerical aperture (NA) of the optical fibre is a measure of light gathering ability of the fibre and is an important parameter which determines the angle of the "light acceptance cone" at the input end of the fibre. The numerical aperture is defined as the sine angle of the acceptance cone. It is given as

$$NA = \sin \theta_a = \sqrt{\mu_f^2 - \mu_c^2}$$

where,  $\mu_f$  and  $\mu_c$  are refractive indices of core and cladding respectively. Thus to determine the numerical aperture, the acceptance angle  $\theta_a$  of the fibre should be known.

In case of experiment setup the numerical aperture can be calculated by the following formulae.

$$\sin \theta = \frac{D/2}{\sqrt{(D/2)^2 + L^2}} = \frac{D}{\sqrt{D^2 + 4L^2}}$$

Calculation

$$\text{Mean value of NA} = \frac{0.8613 + 0.8327 + 0.7682 + 0.7619 + 0.7544}{5}$$

$$NA = \frac{3.9785}{5} = 0.7957$$

Result

NUMERICAL APERTURE OF GIVEN OPTICAL FIBRE:-

$$0.7957$$

where  $D$  is diameter of optical fibre by screen. and  $L$  is distance between optical fibre and screen.

Precautions • Optical fibre should not bend when the light is being propagated.




- Light should be coupled properly into the fibre so that a bright spot of light is obtained.
- Losses due to bending, leakage, dispersion etc should be minimised.
- We should ensure accurate measurement of diameter of the spot.
- Optical fibre axis should be parallel to optical bench.

Sources of error- There may be losses due to bending radiation and leakage, etc.

- The edges on ends of optical fibre might not be exactly sharply cut and the bright spot may not be exactly circular.
- The intensity of light should be measured at the distance between the screen and the output end of the fibre if increased.

Result-

$$\text{Numerical aperture} = 0.7957$$



EXPERIMENT - 6

Aim- To determine the wavelength of sodium light using diffraction-grating

Apparatus- A plane diffraction-grating, spectrometer mercury-comp. prism, spirit-level and magnifying glass.

Theory: The wavelength of any spectral-line can be calculated as:-

$$(a+b) \sin \theta = n\lambda \Rightarrow \lambda = \frac{(a+b) \sin \theta}{n}$$

where  $a+b$  = grating-element,  $\theta$  = angle of diffraction  
 $n$  = order of spectrum

Precautional and sources of errors:-

- ① All adjustments of spectrometer must be correctly done
- ② The readings of the two gratings must be vertical.
- ③ The light must be incident on the unruled-side of the grating. Otherwise refraction occurs at the thick-glass plate after diffraction and the angle measured will not correspond to the diffraction.
- ④ The grating should be as narrow as possible and parallel to the readings.

## Observations & Calculations

Vernier constant of spectrometer =  $1^\circ$

No° of lines per inch of grating  $N = 15000$

Grating element ( $a+b$ ) =  $\frac{25.4}{5000} = 0.00169 \text{ mm.}$

ORDER OF SPECTRUM	VSR Telescope L.H.S	Telescope R.H.S	Half - angular width, $\theta_n$	Average angle, $\theta_n$
0th Order	$180^\circ, 66^\circ$ 364°	$198.8^\circ$ $359.86^\circ$	$U_1$ $U_2$	$\frac{2.58}{2} = 1.29$ $\frac{4.14}{2} = 2.07$
1st Order	$156.66^\circ$ $333.83^\circ$	$1.99^\circ$ $370.66^\circ$	$U_1$ $U_2$	$21.17$ $18.66$

For 1st Order

$n = 1$ ,  $\theta_n = 19.91^\circ$ ,  $\sin \theta_n = 0.340$ ,  $(a+b) = 0.00169 \text{ mm.}$

$$\lambda = (a+b) \sin \theta_n = 0.00169 \times 0.340 \\ = 5755 \text{ Å}$$

Result - Mean value of  $\lambda = 5755 \text{ Å}$

- ⑤ The MFT should be narrow as possible and parallel to the ruling of the grating.
- ⑥ While reading of the telescope, the turntable should be clamped and vice-versa.
- ⑦ Readings should be taken in both the vernier to eliminate the error due to non-coincidence of the centre of the graduated scale with axis of rotation of the spectrometer.
- ⑧ Light should be incident normally on the grating surface. Otherwise the angle of incidence would not be zero which is a condition for the derivation of the formula  $(a+b) \sin \theta = n \lambda$ . Further small change cause a large change in the angle of incidence cause a large error in the angle of diffraction.
- ⑨ Apart from experimental errors, the grating itself may have following errors.
  - (a) Error due to curvature and non-parallelism of the prism. of ~~the~~ the rulings.
  - (b) The error due to variation of the rulings.

## Observation & Calculation:-

Vernier constant of spectrometer =  $1^\circ$ .

No. of lines per inch on grating,  $N = 15000$

$$\text{Grazing element } (a+b) = \frac{25.4}{5000} = 0.00169 \text{ mm}$$

ORDER OF SPECTRUM	USR Telescope L.H.S	Telescope R.H.S		Half angular width $\theta_n$	Average angle $\theta_n$
0th Order	$180.66^\circ$ $364^\circ$	$198.8^\circ$ $359.86^\circ$	$U_1$ $U_2$	$\frac{2.58}{2} = 1.29$ $\frac{4.14}{2} = 2.07$	$1.08$
1st Order	$156.66^\circ$ $333.83^\circ$	$199^\circ$ $370.66^\circ$	$U_1$ $U_2$	$21.17$ $18.66$	$19.91$

For 1st Order:-

$$n=1, \theta_n = 19.91, \sin \theta_n = 0.340, (a+b) = 0.00169 \text{ mm.}$$


$$\lambda = (a+b) \sin \theta_n = 0.00169 \times 0.340.$$

$$= 5755 \text{ Å}$$

Result: Mean value of  $\lambda = 5755 \text{ Å}$ .

(c) The colour of spacing.

Result - Mean value of  $\lambda = 5755 \text{ \AA}$ .



(a) A parallel beam of light falling on the edge of prism (b) measuring angle of prism with spectrometer. (c) Adjustment of the prism on the turn table.

## EXPERIMENT - 7

Aim - To determine the refractive-index of the material of the prism using spectrometer.

Apparatus - A spectrometer, an electric-lamp, a glass prism, a magnifying, an electric lamp and a source of light.

Theory - Let a narrow parallel beam of light coming from the collimator fall on the edge A of prism. It gets partly reflected from it & fall AC. If we place a telescope with its axis parallel to RF or AE and image of the slit screen is seen in it, if we rotate telescope so that its axis becomes parallel to HF or AD, the image of the slit will be seen again in it. The angle of prism A is equal to half angle  $\theta_m$  through which the telescope is rotated i.e.  $\angle A = \frac{\theta_m}{2}$

Formula used - If  $\theta_m$  is the angle of minimum deviation of a ray refracted through a prism, the refractive-index of the material of the prism is

$$\mu = \frac{18n}{\sin \left( \frac{A + \theta_m}{2} \right)}$$

Precautions - ① In finding the angle of the prism, the prism-

## Observation

Value of Vernier scale division  
 Vernier - constant of the spectrometer =  $30''$

S.N.O.	Vernier scale reading M.S.R + (U.S.R X L.C.)		Difference $\theta$ (Degree)	Mean. $\theta_{\text{av}}$ (Degree)	Prism angle $A = \frac{\theta \theta}{2n}$ (degree)
	Telpherite on L.H.S (Right through left + few of prism)	Telpherite on R.H.S (Right through right few of prism)			
(•)	$V_1 = 98.33$ $V_2 = 278.38$ $= 360 - 278.38$ $= 81.62.$	$V_1 = 218.95$ $V_2 = 38.38$	$\theta_1 = 120.12$ $\theta_2 = 120$	$\theta_{\text{av}} = 120.06$ $-$	$A = 60.03.$ $W$

## Angle of minimum Deviation (mm).

Line	Vernier - Scale reading M.S.R + U.S.R X L.C.		Difference $\theta_{\text{av}}$ (degree)	Refractive Index for different colours (n)
	For minimum deviation-ray $V_1$	For direct-ray $V_1$		
Purple	132.5	183.59	51.19	1.650
Indigo	133.375	183.69	50.315	1.6414
Blue	134.35	183.69	49.34	1.6316
Green	135.0833	183.69	48.6067	1.624
Yellow	135.42	183.69	48.27	1.620
Orange	135.5833	183.69	48.1067	1.619
Red	135.75	183.69	47.94.	1.617.

- ② If vernier jaw is at the zero of the scale or telescope is coated, this fact should be taken into account in determining the angle turned.
- ③ Eye piece should be adjusted so that the cross-wire are clearly visible.
- ④ Ensure that the plane of the prism-table is horizontal and axis of the telescope and collimator lie in the same horizontal plane.
- ⑤ The slit of the collimator should be narrow and adjust its position w.r.t the collimator-lens so that the image of the slit is sharply focused.
- ⑥ Both the vernier-scales should be read to avoid errors due to eccentricity.
- ⑦ The telescope should be focused for infinity and collimator is adjusted to give a parallel-beam of light.

#### Sources of error -

- ① Backlash - error may arise during rotation of telescope or prism-table.
- ② Axis of telescope and collimator may not lie in the same horizontal-plane
- ③ The telescope may not be properly set for infinity
- ④ The prism-table may not be exactly horizontal.

Calculation -

Angle of deviation  $\gamma_{\text{violet}} = 51.19$

$$H_{\text{violet}} = \frac{\sin \left( \frac{A + \delta V}{2} \right)}{\sin \left( \frac{A}{2} \right)} = \frac{\sin \left( \frac{60 + 51.19}{2} \right)}{\sin (30^\circ)} = \frac{0.82}{0.50} = 1.64$$

$$\text{Percentage - error} = \frac{1.65 - 1.64}{1.64} \times 100 = 0.6\%$$

Refractive Index (Standard-Value) of the material  
of given-prism from table =  $\frac{H_U + H_R}{2}$ .

$$= \frac{1.650 + 1.617}{2} = 1.6335.$$

## Experiment - 8

Aim - To determine the dispersive-power of the prism using spectrometer and measuring-tube.

Apparatus - Spectrometer, prism, a spirit-level & magnifying-lens, mercury-lamp and sodium-lamp

Theory - The angular dispersion is defined as the product in the extreme-colour of the spectrum.

If  $D_V$  and  $D_R$  are the deviation for violet and red-colour respectively, then angular-dispersion  $\alpha = D_V - D_R$

$$= (\mu_V - 1) A - (\mu_R - 1) A = (\mu_V - \mu_R) A.$$

where  $\mu_V$  and  $\mu_R$  are the refractive indices of violet and red-light,  $A$  is the angle of prism.


The mean-deviation is the deviation produced for the mean-colour of the spectrum. The mean colour is generally taken as the yellow-colour.

$$\text{Mean deviation} = (\mu_Y - 1) A$$

where  $\mu_Y$  is the refractive index of yellow-light. The dispersive power of the material of the prism is defined as the ratio of angular-dispersion to mean deviation.

$$\text{Dispersive-power } w = \frac{(\mu_V - \mu_R) A}{(\mu_Y - 1) A}$$

Aim - To determine the dispersive-power of the prism using spectrometer and mercury source.



$$= \frac{\mu_v - \mu_r}{(\mu_r - 1)}$$

The dispersive-power  $\nu = 1/f$  depends upon the nature of the material material of the prism and is independent of the prism

### Precaution -

- ① In finding the angle of prism, the prism-table should be kept straight
- ② If vernier passes the zero of the scale all the telescope is rotated, then care should be taken in account in determining the angle turned
- ③ Eye-piece should be adjusted so that the lines are clearly visible
- ④ Ensure that the plane of the prism-table is horizontal and the axes of the telescope and collimator lie in same horizontal plane.
- ⑤ The slit of the collimator should be narrow and adjust its position w.r.t the collimator-lens so that the lens image of the slit is sharply focused
- ⑥ Both vernier-scale should be used to avoid error due to eccentricity
- ⑦ The telescope should be focused for infinity and collimator is adjusted to give a parallel beam of light

Observation - Vernier - constant of spectrometer = Angle of prism,  $A$ .

S.NO.	Reading of the vernier V <sub>1</sub>			Reading of vernier V <sub>2</sub>		
	For position of telescope through	Difference (θ)	For position of telescope through	Difference (θ)		
	First-jaw	Second-jaw		First-jaw	Second-jaw	
1	98.33	218.45	120.12	81.07	38.38	120.
2	—	—	—	—	—	—
3	—	—	—	—	—	—

Mean - value of  $\theta = 120.06$ .  
 Angle of prism,  $A = \theta/2 = 60.03$ .

### Angle of minimum - deviation

S.NO.	Line	Vernier-reading	Minimum-deviation ray	Direct Ray	Difference DM
1.	Violet	—	132.5	183.59	Mean DV = 51.19
2.	Yellow	—	135.92	183.69	Mean DV = 48.27
3.	Red	—	135.75	183.69	Mean DR = 47.94

Sources of error:-

- ① Backlash errors may arise during rotations of telescope or prism-table.
- ② Axis of telescope and collimator may not be in the same horizontal plane.
- ③ The telescope may not be properly set for infinity.
- ④ The prism-table itself may not be exactly horizontal.

$$\text{Calculation} - \mu = \frac{n(\mu + Dm)}{n(\mu/2)}$$

$$\mu_0 = \frac{n(\frac{60.03 + 51.19}{2})}{n(60.0312)} = \frac{n(55.61)}{n(30.015)} = \frac{0.8252}{0.5} \\ = 1.6504.$$

$$\mu_r = \frac{n(\frac{60.03 + 48.27}{2})}{n(0.0312)} = 1.62$$

$$\text{Refractive power} = \frac{\mu_0 - \mu_r}{\mu_0 - 1}$$

$$= \frac{1.6504 - 1.616}{1.62 - 1} = \frac{0.0344}{0.02} = 0.05548.$$

Refractive power of material of prism = 0.05548.

## Innovative - Experiment - 1

Aim - To determine the Cauchy's constant of the given prism

Apparatus - Spectrometer prism, Prism clamp, magnifying glass, mercury vapour lamp etc.

Theory - Cauchy's equation is an empirical relationship between the refractive index and wavelength of light for a particular transparent material. The most general form of Cauchy's equation is


$$n(\lambda) = A + \frac{B}{\lambda^2} + \frac{C}{\lambda^4} \quad \text{--- (1)}$$

where  $n$  is refractive index,  $\lambda$  is the wavelength,  $A, B, C$  etc. are coefficient that can be determined for a material by getting the eqn to be measured refractive indices at known-wavelengths.

The refractive index of  $n$  of the material of the prism for a wavelength  $\lambda$  is


$$n = A + \frac{B}{\lambda^2} \quad \text{--- (2)}$$

where  $A$  &  $B$  are called Cauchy's constant for prism. If refractive-index  $n_1$  &  $n_2$  for any two wavelengths  $\lambda_1$  and  $\lambda_2$  are determined by spectrometer, the Cauchy's constant  $A$  and  $B$  can be calculated



Top View of the setup.

Angle of prism




Scale - x-axis : 2 units = 1m<sup>2</sup>  
 y-axis : 2 units = 1 unit  
 of P]

$$\text{Yellow} \rightarrow 2.988 \rightarrow 1/\lambda^2$$

$$\text{Blue} \rightarrow 5.265 \rightarrow 1/\lambda^2.$$

(2.988, 1.620)

(5.265, 1.634)



from the above eq<sup>n</sup>. In particular, the equation is valid for region of normal-dispersion in the visible wavelength region. In the infrared the eq<sup>n</sup> becomes inaccurate, and it can't be representative of regions of anomalous dispersion.

The refractive-index at different wavelength can be calculated as :-

$$\mu = \frac{\sin A_0 + D_m}{2}$$

$$= \sin \frac{A_0}{2}$$

$A_0$  is the angle of prism and  $D_m$  is the angle of minimum deviation of prism.

- Precautionary Precautions:
- ① Care should be taken to ensure proper setting of the spectrometer.
  - ② It should be ensured that the setting of the telescope and collimator aren't disturbed during the course of taking various readings.

#### Source of error -

- ① The slit should be sharp and vertical
- ② The refractive-surfaces of the prism shouldn't be touched with fingers
- ③ The position of angle of minimum-deviation should be accurately determined

Pair of colours	$\lambda_1 \times 10^{-9}$ m	$\lambda_2 \times 10^{-9}$ m	$n_1$	$n_2$	A	B
Yellow & Blue	580	436	1.620	1.6316	1.620	$9.48 \times 10^{-3}$
Green & Violet	546	405	1.624	1.650	1.620	$5.03 \times 10^{-3}$

Calculation - B is the slope of graph.

For B, after plotting graph, according to values,

For 1st pair of light, yellow & blue

$$B = \frac{1.050 - 1.024}{6.103 - 3.353} = \frac{0.026}{2.75} = 9.43 \times 10^{-3} = B$$

$n = A + \frac{B}{\lambda^2}$ ,  $A = 1.620$  or  $\frac{B}{\lambda^2}$  is too small

For 2nd pair of light, green & violet

$$B = \frac{1.6316 - 1.620}{5.265 - 2.988} = \frac{0.0116}{2.277} \\ = 5.09 \times 10^{-3}$$

$n = A + \frac{B}{\lambda^2}$ ,  $A = 1.620$  or  $\frac{B}{\lambda^2}$  is too small.

↑

Scale: x-axis : 2 units =  $1\text{m}^2$   
 y-axis : 2 units =  $10\text{nt}$   
 or  $0.1$

Green  $\rightarrow 3.353 \rightarrow 11\lambda^2$

Violet  $\rightarrow 6.103 \rightarrow 11\lambda^2$

n1 3


2

3.353, 1.624

R 6.103, 1.650

1

u



## Innovative Experiment - 2

Aim - To determine the dispersive power of a plane transmission grating arranged for normal incidence.

Apparatus - spectrometer, given - grating, mercury vapour-lamp etc.

Principle - The dispersive power of grating, is the ratio of change in angle of diffraction to the corresponding change in wavelength of any two neighbouring rays. Let  $\theta$  wavelengths  $\lambda$  and  $d\lambda$  be diffracted through  $\Theta$  and  $d\Theta$ .

$$\text{Dispersive power of grating} = d\Theta/d\lambda$$

For a grating for normal incidence,

$$\sin \Theta = N \sin \Theta \quad \text{--- (1)}$$

where  $\Theta$  is the angle of diffraction,  $N$  is order of spectrum.  $N$  is the no. of lines per unit length and  $\lambda$  is the wavelength of light.


Differentiating, we get  $-N d\Theta = N d\lambda$ .

$$\frac{d\Theta}{d\lambda} = \frac{N}{\cos \Theta}$$

Result -

(1) Mean dispersive power of grating for violet region = 636.44 degrees/m

Aim - To find the diffraction for a grating.



② Mean dispersive-power of grating in yellow - region  
= 636.94 degree lm.

### Precautions -

- ① Care should be taken to ensure proper setting of the spectrometer and these settings of the telescope and the collimator aren't touched during the course of taking other various readings.
- ② The position of the grating adjusted to be normal to the incoming-light from collimator should n't be disturbed.
- ③ It is necessary to point the slit toward the brightest - part of the source in order to obtain measurable intensity of the lines of different - colour especially in high order spectra which sharply with decrease in order.

## Observation & Calculations

Vernier constant of spectrometer = 10

No. of lines per inch on grating,  $N = 1500$ .

$$\text{Grating element (a+b)} = \frac{25.4}{5000} = 0.00169 \text{ mm}$$

ORDER OF SPECTRUM	VSR Tellerator LHS	Tellerator RHS		Half angular width, $\theta_n$	Average angle, $\theta_n$
0 <sup>th</sup> Order	180.66° 364°	198.8° 359.86°	v <sub>1</sub> v <sub>2</sub>	$\frac{2.58}{2} = 1.29$ $\frac{4.14}{2} = 2.07$	1.08.
1 <sup>st</sup> Order	156.66° 333.83°	199° 370.66°	v <sub>1</sub> v <sub>2</sub>	21.17 18.66	19.91

Calculation (a+b)  $\sin \theta = n\lambda$

$$(a+b) \cos \theta = n \frac{d\lambda}{d\theta} \Rightarrow \frac{d\theta}{d\lambda} = \frac{n}{(a+b) \cos \theta}$$

$$\text{For } n=0, \frac{d\theta}{d\lambda} = 0$$

$$\begin{aligned} \text{For } n=1, \frac{d\theta}{d\lambda} &= \frac{1}{(a+b) \cos \theta} \\ &= \frac{1}{(a+b) \cos (20.8)} \\ &= \frac{1}{0.00169 \times 0.934825} \\ &= \frac{1}{0.00157} = 636.942675159 \end{aligned}$$

Result - Dispersion power of grating is  
636.94 degree/m.