1. The percentage cum position of mixture of NaOH and NaCl, 10gm of which has been dissolved litre (given N/10 HCl)

Titration Type	Indicator Name	Properties of Indicator	Significance of Titration	Reaction Involved	Solution in Burette	Solution in Conical Flask	Observations	Extra Info
Acid- Base	Methyl Orange	pH range: 3.1-4.4	Indicates endpoint by changing color from red to yellow	HCI + NaOH → NaCI + H2O	N/10 HCI solution	NaOH and NaCI mixture solution	The initial color of the solution is red, which changes to yellow upon addition of the titrant. The titration is complete when the color changes from red to yellow.	

2. To find the amount of oxalic acid dissolved per litre. You are provided with N/10 NaOH solution

Titration Type	Indicator Name	Properties of Indicator	Significance of Titration	Reaction Involved	Solution in Burette	Solution in Conical Flask	Observations	Extra Info
Acid- Base	Phenolphthalein	pH range: 8.2-10.0	Indicates endpoint by changing color from colorless to pink	NaOH + H2C2O4 → Na2C2O4 + 2H2O	N/10 NaOH solution	Given oxalic acid solution	The initial color of the solution is colorless, which changes to pink upon addition of the titrant. The titration is complete when the color changes from colorless to pink.	-

3. To find amount the amount of oxalic acid dissolved per litre, You are provided with N/10 KMnO4 solution (using self-indicator)

Titration Type	Indicator Name	Properties of Indicator	Significance of Titration	Reaction Involved	Solution in Burette	Solution in Conical Flask	Observations	Extra Info
Oxidation- Reduction	Self- Indicator (Oxalic acid)	A weak acid that changes color as the pH changes	Indicates endpoint by the color change from pink to colorless	H2C2O4 + KMnO4 → CO2 + MnO2 + H2O	N/10 KMnO4 solution	Given oxalic acid solution	The initial color of the solution is pink, which fades as the titrant is added. At the endpoint, the color changes from pink to colorless.	Self- indicator used in this titration is oxalic acid itself

HERE AT LAST LIGHT PINK COLOUR IS OBTAINED INSTEAD OF COLOURLESS. PLEASE NOTE in 3)

4. MOHR'S METHOD

To determine the strength of chloride ion, being provided with approx N/50 silver nitrate Solution (Mohr's METHOD)

Titration Type	Indicator Name	Properties of Indicator	Significance of Titration	Reaction Involved	Solution in Burette	Solution in Conical Flask	Extra Info
Precipitation titration	Potassium chromate (K2CrO4)	Yellow in acidic medium and red in alkaline medium	Appearance of red precipitate indicates the endpoint	AgNO3 + NaCl → AgCl + NaNO3	Approx. N/50 AgNO3 solution	Chloride ion containing solution	Potassium chromate acts as an indicator for the formation of silver chromate, which is used as an endpoint indicator. The endpoint is marked by the appearance of a red precipitate of silver chromate.

5. To determine the strength of given CuSO4 solution, provided with N/50 sodium thiosulphate solution

Titration Type	Indicator Name	Properties of Indicator	Significance of Titration	Reaction Involved	Solution in Burette	Solution in Conical Flask	Extra Info
Redox	Starch	Sensitive to oxidation and reduction reactions	End point: Disappearance of blue colour	CuSO4 + Na2S2O3 → CuS2O3 + Na2SO4	N/50 Na2S2O3	CuSO4 solution	The solution in the conical flask must be acidified with dilute H2SO4 to prevent the formation of Cu(OH)2

KI IS ADDED IN THIS SOLUTION SUCH THAT COPPER SULPHATE CONVERTS TO CU212 (WHITE PPT)

AND EQUIVALENT AMOUNT OF IODINE IS FORMED

7. To determine the total hardness of given water sample by EDTA method.

Titration Type	Indicator Name	Properties of Indicator	Significance of Titration	Reaction Involved	Solution in Burette	Solution in Conical Flask	Extra Info
Complexometric Titration	Eriochrome Black T	Wine red in presence of Ca and Mg, blue in absence	Endpoint detection	Ca2+ + Mg2+ + EDTA -> CaEDTA + MgEDTA	EDTA Solution	Water Sample	Adjust the pH of the water sample to around 10 using ammonia buffer.

EBT IS WEAK LIGAND AND COMPLEX WITH EBT IS COLOURED

*REVISE STRUCTURES

$$\begin{array}{l} {{\text{HOOC-CH}_{_2}}} \\ {{\text{HOOC-CH}_{_2}}} > {N{\text{-CH=CH-N}}} < \frac{{{\text{CH}_{_2}}{\text{-COOH}}}}{{{\text{CH}_{_2}}{\text{-COOH}}}} \end{array}$$

EDTA

$$\begin{array}{c|c}
& & & \\
& & & \\
& & & \\
& & & \\
& & & \\
& & & \\
& & & \\
& & & \\
& & & \\
& & & \\
& & & \\
& & & \\
& & & \\
& & & \\
& & & \\
& & & \\
& & & \\
& & & \\
& & & \\
& & & \\
& & & \\
& & & \\
& & & \\
& & & \\
& & & \\
& & & \\
& & & \\
& & & \\
& & & \\
& & & \\
& & & \\
& & & \\
& & & \\
& & & \\
& & & \\
& & & \\
& & & \\
& & & \\
& & & \\
& & & \\
& & & \\
& & & \\
& & & \\
& & & \\
& & & \\
& & & \\
& & & \\
& & & \\
& & & \\
& & & \\
& & & \\
& & & \\
& & & \\
& & & \\
& & & \\
& & & \\
& & & \\
& & & \\
& & & \\
& & & \\
& & & \\
& & & \\
& & & \\
& & & \\
& & & \\
& & & \\
& & & \\
& & & \\
& & & \\
& & & \\
& & & \\
& & & \\
& & & \\
& & & \\
& & & \\
& & & \\
& & & \\
& & & \\
& & & \\
& & & \\
& & & \\
& & & \\
& & & \\
& & & \\
& & & \\
& & & \\
& & & \\
& & & \\
& & & \\
& & & \\
& & & \\
& & & \\
& & & \\
& & & \\
& & & \\
& & & \\
& & & \\
& & & \\
& & & \\
& & & \\
& & & \\
& & & \\
& & & \\
& & & \\
& & & \\
& & & \\
& & & \\
& & & \\
& & & \\
& & & \\
& & & \\
& & & \\
& & & \\
& & & \\
& & & \\
& & & \\
& & & \\
& & & \\
& & & \\
& & & \\
& & & \\
& & & \\
& & & \\
& & & \\
& & & \\
& & & \\
& & & \\
& & & \\
& & & \\
& & & \\
& & & \\
& & & \\
& & & \\
& & & \\
& & & \\
& & & \\
& & & \\
& & & \\
& & & \\
& & & \\
& & & \\
& & & \\
& & & \\
& & & \\
& & & \\
& & & \\
& & & \\
& & & \\
& & & \\
& & & \\
& & & \\
& & & \\
& & & \\
& & & \\
& & & \\
& & & \\
& & & \\
& & & \\
& & & \\
& & & \\
& & & \\
& & & \\
& & & \\
& & & \\
& & & \\
& & & \\
& & & \\
& & & \\
& & & \\
& & & \\
& & & \\
& & & \\
& & & \\
& & & \\
& & & \\
& & & \\
& & & \\
& & & \\
& & & \\
& & & \\
& & & \\
& & & \\
& & & \\
& & & \\
& & & \\
& & & \\
& & & \\
& & & \\
& & & \\
& & & \\
& & & \\
& & & \\
& & & \\
& & & \\
& & & \\
& & & \\
& & & \\
& & & \\
& & & \\
& & & \\
& & & \\
& & & \\
& & & \\
& & & \\
& & & \\
& & & \\
& & & \\
& & & \\
& & & \\
& & & \\
& & & \\
& & & \\
& & & \\
& & & \\
& & & \\
& & & \\
& & & \\
& & & \\
& & & \\
& & & \\
& & & \\
& & & \\
& & & \\
& & & \\
& & & \\
& & & \\
& & & \\
& & & \\
& & & \\
& & & \\
& & & \\
& & & \\
& & & \\
& & & \\
& & & \\
& & & \\
& & & \\
& & & \\
& & & \\
& & & \\
& & & \\
& & & \\
& & & \\
& & & \\
& & & \\
& & & \\
& & & \\
& & & \\
& & & \\
& & & \\
& & & \\
& & & \\
& & & \\
& & & \\
& & &$$

EBT STRUCTURE

7. To determine the strength of magnesium ions in given MgSO4 solution by complexometric titration

Titration Type	Indicator Name	Properties of Indicator	Significance of Titration	Reaction Involved	Solution in Burette	Solution in Conical Flask	Extra Info
Complexometric	Eriochrome Black T	Forms a wine-red complex with metal ions	End point: Colour change from wine- red to blue	Mg2+ + EDTA → MgEDTA2-	N/50 EDTA	MgSO4 solution + NH4Cl buffer	The pH of the solution should be maintained between 10 and 11 for accurate results

EXPERIMENT 8

Here is the corrected table for the titration to determine the alkalinity of a water sample using phenolphthalein and methyl orange:

Titration Type	Indicator Name	Properties of Indicator	Significance of Titration	Reaction Involved	Solution in Burette	Solution in Conical Flask	Extra Info
Alkalinity determination	Phenolphthalein	Colorless in acidic pH, pink in basic pH	Marks the endpoint of the first part of the titration when the sample changes from acidic to basic	H2CO3 + 2NaOH -> Na2CO3 + 2H2O	N/10 HCI	Water Sample	-
Alkalinity determination	Methyl Orange	Yellow in acidic pH, red in basic pH	Marks the endpoint of the second part of the titration when the sample changes from basic to alkaline	HCO3- + OH > CO3-2 + H2O	N/10 HCI	Water Sample + Phenolphthalein	-

VIVA QUESTIONS

1. Name All types of titration.

- 1. NaOH and NaCl mixture: Acid-base titration.
- 2. Determination of oxalic acid using NaOH: Acid-base titration
- 3. Determination of oxalic acid using KMnO4: Redox titration
- 4. Determination of chloride ion using silver nitrate (Mohr's Method): Precipitation titration
- 5. Determination of CuSO4 using sodium thiosulfate: Redox titration
- 6. Determination of total hardness of water using EDTA: Complexometric titration
- 7. Determination of Ca2+ using EDTA: Complexometric titration
- 8. Determination of Mg2+ using EDTA: Complexometric titration
- 9. Determination of alkalinity using phenolphthalein and methyl orange: Acidbase titration

2. What is autocatalyst?

An autocatalyst is a substance that catalyzes a chemical reaction and is itself one of the reactants or products of the reaction. In other words, the substance acts as a catalyst to speed up the reaction while being consumed or produced in the process. As the reaction proceeds, the concentration of the autocatalyst increases, which in turn accelerates the reaction even further. This positive feedback loop can lead to a rapid increase in the reaction rate, and in some cases, it can lead to a runaway reaction or an explosion. Autocatalysis is a common phenomenon in many chemical and biochemical systems, including enzyme-catalyzed reactions, polymerization reactions, and oscillating reactions.

3 Give the structure of EBT

EBT stands for Eriochrome Black T and it is a complexometric indicator used in analytical chemistry

The sulfonate and carboxylic acid groups in EBT are responsible for complexing with metal ions, and the heterocyclic ring imparts color to the compound. In its protonated form, EBT is yellow, and when it forms a complex with metal ions, it turns wine-red.

$$\begin{array}{c|c} & & & & \\ & & & & \\ & & & \\ & & & \\ & & & \\ & & & \\ & \\ & & \\ & & \\ & & \\ & & \\ & & \\ & & \\ & & \\ & & \\ & & \\ & & \\ & & \\ & & \\ & & \\ &$$

EBT STRUCTURE

What is hardness?

Hardness in water refers to the amount of dissolved minerals, primarily calcium and magnesium, present in the water.

These minerals can make the water "hard" and can cause problems such as scale buildup in pipes and appliances, soap scum, and reduced lathering of soaps and detergents. Hardness is typically measured in terms of the amount of calcium carbonate present in the water, and is often expressed in units of milligrams per liter (mg/L) or parts per million (ppm). Hardness can be classified as either temporary or permanent, depending on the type of minerals present in the water. Temporary hardness can be removed by boiling the water, while permanent hardness requires other treatment methods.

What is a buffer solution?

A buffer solution is a solution that can resist changes in pH when small amounts of an acid or a base are added to it. It typically consists of a weak acid and its corresponding conjugate base, or a weak base and its corresponding conjugate acid. Buffers are commonly used in laboratory experiments to maintain a stable pH environment for chemical reactions. For example, biological systems rely on buffers to maintain a stable pH environment for proper functioning.

Can NaCl produce hardness on dissolving it in water?

No, NaCl (sodium chloride) does not produce hardness when dissolved in water. Hardness in water is mainly caused by the presence of dissolved calcium and magnesium ions. NaCl is a salt that dissociates in water to form sodium ions and chloride ions, which do not contribute to water hardness.

NAME THE SALTS WHICH FORM PERMANENT HARDNESS?

Some of the salts that can form permanent hardness in water are:

- 1. Calcium sulfate (CaSO4)
- 2. Calcium chloride (CaCl2)
- 3. Magnesium sulfate (MgSO4)
- 4. Magnesium chloride (MgCl2)
- 5. Iron(III) sulfate (Fe2(SO4)3)

GIVE definition of all TITRATION?

- 1. Acid-Base Titration: This experiment involves the titration of a strong acid (e.g. hydrochloric acid) with a strong base (e.g. sodium hydroxide) or vice versa. The aim is to determine the concentration of the acid or base by measuring the amount of titrant required to neutralize the solution.
- 2. Iodometric Titration: This experiment involves the titration of iodine with a reducing agent such as thiosulfate. The aim is to determine the concentration

- of the oxidizing agent by measuring the amount of reducing agent required to react with the iodine.
- 3. Redox Titration: This experiment involves the titration of a reducing agent with an oxidizing agent or vice versa. The aim is to determine the concentration of the oxidizing or reducing agent by measuring the amount of titrant required to react with the analyte.
- 4. Precipitation Titration: This experiment involves the titration of an analyte with a titrant that causes a precipitate to form. The aim is to determine the concentration of the analyte by measuring the amount of titrant required to cause the precipitation.
- 5. Complexometric Titration: This experiment involves the titration of an analyte with a chelating agent such as EDTA. The aim is to determine the concentration of the analyte by measuring the amount of chelating agent required to form a stable complex with the analyte.

What is the difference between end point and equivalence point?

The endpoint is the point in a titration where a chemical indicator changes color, signaling the completion of the reaction. It is a visual indication of the equivalence point, which is the theoretical point in the titration where the number of moles of the reactant being titrated is equal to the number of moles of the titrant added. While the endpoint is observed through the use of an indicator, the equivalence point is calculated based on the stoichiometry of the reaction. The difference between the endpoint and equivalence point is that the endpoint is a physical observation while the equivalence point is a theoretical calculation.

Suppose four beakers are given one contains calcium carbonate, the second one contains calcium sulphate, the third one contains calcium hydrogen carbonate and the fourth contains iron. Identify the type of hardness present in each.

Based on the given information, we can identify the type of hardness present in each beaker as follows:

- 1. Beaker containing calcium carbonate: This type of hardness is known as temporary hardness, which is caused by the presence of dissolved bicarbonate and carbonate ions in water. Calcium carbonate is one of the minerals that can contribute to temporary hardness.
- 2. Beaker containing calcium sulphate: This type of hardness is known as permanent hardness, which is caused by the presence of dissolved sulfate and chloride ions in water. Calcium sulphate is one of the minerals that can contribute to permanent hardness.
- 3. Beaker containing calcium hydrogen carbonate: This type of hardness is also temporary hardness, which is caused by the presence of dissolved bicarbonate

- and carbonate ions in water. Calcium hydrogen carbonate is a form of temporary hardness that is present in some natural waters.
- 4. Beaker containing iron: Iron does not contribute to water hardness.

Is edta soluble in water? Why?

Yes, EDTA (ethylenediaminetetraacetic acid) is soluble in water. This is because EDTA is a polar molecule, meaning it has a positive and negative end, which allows it to interact with water molecules through hydrogen bonding. Additionally, EDTA is a weak acid and forms negatively charged ions in solution, which also contributes to its solubility in water.

It is a water-soluble organic compound that can form complexes with metal ions, making it a commonly used chelating agent in various applications, including analytical chemistry and water treatment.

What is alkalinity?

Alkalinity refers to the capacity of water to resist a change in pH when an acid is added to it. It is primarily due to the presence of bicarbonate, carbonate, and hydroxide ions in the water. Alkalinity is an important water quality parameter, especially in water treatment processes and aquatic ecosystems.

Examples of alkaline substances that can contribute to alkalinity include bicarbonates in baking soda and carbonates in limestone. In natural water sources, such as rivers and lakes, dissolved minerals like calcium, magnesium, and bicarbonates contribute to alkalinity. Alkalinity can also be artificially increased in water treatment processes by adding chemicals like lime or sodium carbonate.

What is dimineralization?

Deminerlization (also spelled demineralization) is the process of removing minerals, salts, and other impurities from water. It involves the use of special resins or ion-exchange membranes that selectively remove ions and minerals from the water, producing purified water that is free from dissolved solids. Demineralization is commonly used in industries such as pharmaceuticals, electronics, and power generation, where the presence of impurities can be harmful to processes and equipment.

What is sodium thiosulphate and indicator used in iodometric titration?

In iodometric titration, sodium thiosulfate (Na2S2O3) is used as a reducing agent to convert iodine (I2) to iodide ions (I-) in the presence of an acid. The indicator used in iodometric titration is starch, which forms a blue-black complex with iodine.

What is permanent hardness and temporary hardness and explain ways to remove them?

Hardness in water is due to the presence of calcium and magnesium ions, which can form insoluble salts with soap and cause scaling in pipes and water heaters. Hardness is classified into two types: permanent hardness and temporary hardness.

Temporary hardness is caused by the presence of bicarbonate ions (HCO3-) in water, which react with calcium and magnesium ions to form insoluble carbonates that can be removed by boiling the water or adding lime (calcium hydroxide) to precipitate the carbonates. This type of hardness can also be removed by using an ion-exchange resin.

Permanent hardness is caused by the presence of other soluble calcium and magnesium salts such as sulfates, chlorides, and nitrates, which cannot be removed by boiling. This type of hardness can be removed by using a water softener, which contains a resin that exchanges sodium ions for the calcium and magnesium ions.

In summary, temporary hardness can be removed by boiling or adding lime, while permanent hardness requires the use of a water softener or an ion-exchange resin.

In which form EDTA is used and why?

EDTA (Ethylenediaminetetraacetic acid) is used in the form of its sodium salt, known as disodium EDTA or Na2EDTA. This is because disodium EDTA is highly watersoluble, which makes it easier to use in aqueous solutions. In addition, disodium EDTA has a higher stability constant than EDTA itself, which means it can form more stable complexes with metal ions.

EDTA is a chelating agent, which means it can form complexes with metal ions. In analytical chemistry, EDTA is used to determine the concentration of metal ions in a sample through complexometric titrations. In addition, EDTA is used in various industrial applications, such as in water treatment to remove metal ions and as a stabilizer in cosmetics and food products.

How is brackish water desalinised?

Brackish water, which has a higher salinity level than freshwater but lower salinity than seawater, can be desalinated using various methods such as reverse osmosis, nanofiltration, electrodialysis, and distillation.

One common method for desalinating brackish water is reverse osmosis. In this process, brackish water is pressurized and passed through a semi-permeable membrane that allows water molecules to pass through but blocks salt and other impurities. The purified water is then collected on the other side of the membrane. The remaining brine solution containing the removed salts and impurities is then discharged or further treated for safe disposal.

Another method is electrodialysis, which uses an electric field to separate ions from water molecules. In this process, the brackish water is passed through a stack of alternating positively and negatively charged membranes. The charged membranes attract and remove the ions, leaving behind desalinated water.

Distillation is also a common method used for desalinating brackish water. In this process, the brackish water is boiled and the resulting steam is condensed to produce fresh water. The salts and other impurities are left behind in the boiling chamber.

Each of these methods has its advantages and disadvantages and is selected based on factors such as cost, energy requirements, and the quality of the treated water needed.

WHY TITRATION PROCESS IS INACCURATE/ERROR CAUSES?

Titration can be inaccurate due to various reasons:

- 1. Human error: The person performing the titration may make errors in reading the burette, measuring the volume of the solution, or mixing the reagents.
- 2. Impure chemicals: If the chemicals used in the titration are not pure, it can affect the accuracy of the results.
- 3. Air bubbles: Air bubbles in the burette, pipette or conical flask can affect the volume of the solution, leading to inaccurate results.
- 4. Contamination: If the solution being titrated is contaminated with other substances, it can affect the reaction and lead to inaccurate results.
- 5. Improper calibration: If the equipment used in the titration is not calibrated properly, it can result in inaccurate readings.
- 6. End-point determination: Determining the end-point of the titration can be subjective and can differ from person to person.

- 7. Temperature changes: Changes in temperature can affect the volume and concentration of the solution, leading to inaccurate results.
- 1. Molarity (M): Molarity is defined as the number of moles of solute per liter of solution.

Molarity = (moles of solute) / (volume of solution in liters)

Unit: mol/L

1. Normality (N): Normality is defined as the number of equivalents of solute per liter of solution. It is used for acid-base and redox titrations.

Normality = (number of equivalents of solute) / (volume of solution in liters)

Unit: eq/L

2. Molality (m): Molality is defined as the number of moles of solute per kilogram of solvent.

Molality = (moles of solute) / (mass of solvent in kg)

Unit: mol/kg

- 3. Percentage (%) solutions: Percentage solutions are expressed as the mass of solute per 100 parts of solution.
- % Concentration = (mass of solute / mass of solution) x 100%
 - 4. Calcium Carbonate Equivalent (CCE): Calcium Carbonate Equivalent is a measure of the amount of acid-neutralizing capacity of an unknown material. It is used to express the total concentration of hardness ions in water.

CCE = (mg of hardness ion) / (mg of CaCO3)

Unit: mg/L or ppm (parts per million)

which indicator is used for determining alkalinity?

The commonly used indicators for determining alkalinity are phenolphthalein and methyl orange.

What is the use of H2SO4 in formation of Urea Formaldehyde?

Role 1 :- H2SO4 is used as a catalyst for ionization of formaldehyde in methalol carbonation.

Role 2:- H2SO4 provide the he acidic conditions required for the synthesis of UF resin

Role 3:- H2SO4 is a very good dehydrating agent, so it absorbs water and prevent the further dilution of the reactant (formaline and urea)

Use of h2so4 in case of chemical reactions?

Sulfuric acid is a strong acid and has many uses in various chemical reactions. Some of its common uses are:

- 1. Esterification: Sulphuric acid is used to catalyze the reaction between carboxylic acids and alcohols, forming esters.
- 2. Dehydration: Sulphuric acid can act as a dehydrating agent by removing water molecules from organic compounds, such as alcohols and sugars.
- 3. Polymerization: In polymerization reactions, sulphuric acid is used to catalyze the reaction between monomers to form polymers.
- 4. Neutralization: Sulphuric acid is often used to neutralize basic substances and adjust pH levels in various chemical processes.
- 5. Drying agent: It is also used as a drying agent in laboratories to remove water from organic solvents.
- 6. Electrolyte: Sulphuric acid is used as an electrolyte in lead-acid batteries, which are commonly used in vehicles and as backup power supplies.

what is the difference between primary and secondary standard in titration?

In titration, a primary standard is a substance that is pure, stable, and has a high molecular weight. It can be accurately weighed and dissolved in water to prepare a standard solution of known concentration that can be used to determine the concentration of an unknown solution through titration.

On the other hand, a secondary standard is a substance that is not as pure as a primary standard and has a lower molecular weight. It is not used to prepare a standard solution directly but is instead standardized against a primary standard. In other words, the concentration of a secondary standard is determined by comparing its reaction with a known quantity of a primary standard.

The use of primary standards ensures high accuracy in titration while secondary standards may have a degree of uncertainty due to impurities.

why do we heat the solution after using h2so4?

Heating the solution after adding sulfuric acid (H2SO4) serves several purposes in some chemical reactions, including:

- 1. To remove any excess water: In some reactions, adding H2SO4 to the solution can cause water to be produced. By heating the solution, any excess water can evaporate, leaving only the desired products.
- 2. To speed up the reaction: Heating can increase the rate of many chemical reactions, so by heating the solution, the reaction can proceed more quickly.
- 3. To decompose some compounds: Some compounds may be unstable at room temperature and need to be heated to decompose them into more stable products.
- 4. To dissolve a solid: If a solid reactant is added to the solution, heating can help dissolve it more quickly and thoroughly, ensuring that it reacts completely with the other reagents.

Overall, heating the solution after adding H2SO4 can help ensure that the reaction proceeds as desired and that the desired products are formed.

why do we use starch in iodometric titration?

Starch is used as an indicator in iodometric titration to detect the endpoint of the reaction between iodine and thiosulfate. When iodine is added to the solution, it reacts with the thiosulfate to form a colorless complex. As the concentration of iodine decreases, the complex breaks down and iodine begins to form again. At the endpoint, all of the iodine has reacted with the thiosulfate, and the complex breaks down completely, resulting in the formation of a blue-black starch-iodine complex. The appearance of this blue-black color indicates that all of the iodine has been consumed, and the endpoint of the titration has been reached.

what salt of EDTA DO WE USE SINCE EDTA IS INSOLUBLE IN WATER?

EDTA (ethylene-diamine-tetra-acetic acid) itself is insoluble in water. Therefore, its disodium salt, Na2EDTA, is typically used in complexometric titrations as it is highly soluble in water.

Chem Viva Questions

what is phenolpthalein colour

Phenolphthalein is a colorless chemical compound, which turns pink or fuchsia in the presence of an alkaline solution with a pH greater than 8.2. It is often used as an indicator in titrations to determine the endpoint of a reaction.

what is N=N in Erichrome Black T?

Eriochrome Black T (EBT) is an indicator used in complexometry. It contains an N=N functional group, which is the nitrogen-nitrogen double bond. The use of Eriochrome Black T in complexometry is based on the fact that it forms highly colored chelates with metal ions, which makes it useful for the determination of the presence and concentration of metal ions in solution.

what is h2so4 use in formation of urea formaldehyde?

In the synthesis of Urea-Formaldehyde (UF) resins, H2SO4, also known as sulfuric acid, is used as a catalyst. Sulfuric acid helps to speed up the reaction between urea and formaldehyde to form the UF resin. The acid catalyst acts by protonating the nitrogen atom in the urea molecule, making it more reactive and increasing the rate of reaction with formaldehyde. The reaction between urea and formaldehyde in the presence of sulfuric acid forms a condensation product, which polymerizes to form the UF resin.

what is the difference between end point and equivalence point in titration?

In titration, the endpoint and equivalence point refer to different stages of the reaction between the titrant (the solution being added to the analyte) and the analyte (the solution being titrated).

The endpoint is the stage of the titration at which a visible change in the solution occurs, indicating that the reaction is complete. This is often marked by a color change in an indicator. The endpoint is used to approximate the equivalence point.

The equivalence point is the stage of the titration at which all the reactants have reacted and the reaction has gone to completion. It is the point at which the number of moles of titrant added is equal to the number of moles of analyte present in the solution. The equivalence point is often determined by measuring the change in pH, conductivity, or potential of the solution during the titration.

In summary, the endpoint is a visible indication of the progress of the reaction, while the equivalence point is the exact point at which the reaction is complete.

Which indicator is used to determine alkalinity of water due to OH- phenolphthalein or Methyl Orange?

Methyl Orange is commonly used to determine the alkalinity of water, while Phenolphthalein is commonly used to determine the acidic or basic nature of a solution.

Methyl Orange changes color from yellow to red as the pH of the solution increases from acidic to slightly basic (pH 3.1 to 4.4). This makes it a useful indicator for determining the presence and concentration of hydroxide ions (OH-) in a solution, which are a measure of its alkalinity.

Phenolphthalein changes color from colorless to pink or fuchsia as the pH of the solution increases from acidic to basic (pH 8.2 and above). This makes it a useful indicator for determining the presence and concentration of hydronium ions (H3O+) in a solution, which are a measure of its acidity.

In conclusion, Methyl Orange is the preferred indicator for determining the alkalinity of water due to hydroxide ions, while Phenolphthalein is the preferred indicator for determining the acidic or basic nature of a solution.

why do we make solutions in distilled water?

Distilled water is commonly used to make solutions because it is pure and free of impurities that may interfere with the intended reaction or experiment.

Distilled water is produced by boiling tap water and condensing the steam, which leaves behind dissolved minerals and other impurities. This makes distilled water an ideal choice for preparing solutions that require a low concentration of dissolved solids or a specific pH.

For example, in chemical experiments or analytical procedures, using distilled water can prevent contamination of the solution and ensure accurate results. In biological experiments, distilled water is used to make solutions for cell cultures or tissue culture media because it does not contain the minerals or organic compounds that may be toxic to the cells.

In conclusion, distilled water is used to make solutions because it is pure and free of impurities, making it ideal for use in a variety of applications where a low concentration of dissolved solids or a specific pH is required.

why is h2so4 used in various reactions?

Sulfuric acid (H2SO4) is commonly used in titrations as a catalyst because it accelerates the reaction between the titrant (the solution being added to the analyte) and the analyte (the solution being titrated). By increasing the rate of reaction, sulfuric acid makes it easier to detect the endpoint of the titration, when the reaction has gone to completion.

In addition, sulfuric acid can also help to maintain the pH of the solution at a desired level during the titration. This is important because the rate of the reaction between the titrant and analyte is often pH-dependent. By maintaining a constant pH, sulfuric acid helps to ensure that the reaction occurs at a predictable rate, making it easier to determine the endpoint of the titration.

In conclusion, sulfuric acid is used in titrations as a catalyst because it accelerates the reaction between the titrant and analyte and helps to maintain the pH of the solution during the titration, making it easier to determine the endpoint of the reaction.

why distilled water is used for rinsening?

Distilled water is commonly used for rinsing because it is pure and free of impurities that may interfere with the intended use or experiment.

Rinsing with distilled water can help to remove any residual impurities or contaminants that may be present on a surface or in a solution. This is important in many applications, such as rinsing laboratory glassware or cleaning delicate electronic components, where residual impurities can cause contamination, corrosion, or other undesirable effects.

In addition, distilled water is also ideal for rinsing because it does not contain dissolved minerals or other impurities that may leave deposits on the surface being rinsed. This is particularly important in applications where the residue could interfere with subsequent reactions or experiments.

In conclusion, distilled water is used for rinsing because it is pure and free of impurities, making it ideal for removing residual contaminants and preventing deposits from forming on the surface being rinsed.

what is the indicator used in cuso4 titration?

In the titration of copper(II) sulfate (CuSO4) with a base, a suitable indicator must be used to determine the endpoint of the reaction. There are several indicators that can be used for this purpose, including:

- 1. Phenolphthalein: Phenolphthalein is a commonly used indicator for titrations involving strong acids and strong bases. In a neutral or acidic solution, phenolphthalein is colorless. As the solution becomes basic, the color changes from colorless to pink.
- 2. Methyl Orange: Methyl orange is a pH indicator that is commonly used in titrations involving weak acids and strong bases. In acidic solutions, methyl orange is yellow, and in basic solutions, it is orange.
- 3. Bromocresol Green: Bromocresol green is a pH indicator that is also commonly used in titrations involving weak acids and strong bases. In acidic solutions, bromocresol green is yellow, and in basic solutions, it is blue.

In conclusion, phenolphthalein, methyl orange, and bromocresol green are some of the indicators that can be used in the titration of copper(II) sulfate with a base. The choice of indicator will depend on the specific requirements of the titration, such as the desired endpoint and the sensitivity of the indicator.

What is Auto Catalysis?

Autocatalysis is a chemical reaction in which one of the reaction products acts as a catalyst, accelerating its own formation in a positive feedback loop. This can lead to a rapid increase in the concentration of the reaction product, potentially resulting in a self-sustaining reaction.

What is Self Indicator?

A self-indicator is a substance that can indicate its own concentration or presence in a solution without the need for an external reagent or reference. It is a type of chemical indicator that changes its physical or chemical properties in response to the concentration of the substance being measured. For example, a pH indicator can change color in response to the pH of the solution it is in, allowing the pH to be determined by observing the color change. Self-indicators can be useful in a variety of applications, including monitoring the progress of chemical reactions, determining the presence of specific substances in a sample, and quality control in industries such as food and pharmaceuticals.

What is the Unit of Strength?

Strength of a solution:

• "The strength of a solution is defined as the mass of solute (in grams) dissolved in one liter of solution."

Formula: $Strength = \frac{Massof soluteing rams}{Volume of solution in litre}$

- It can be expressed as:
- The units of strength are gL-1 or gdm-3.

Give an example of precipitation titration. Also give its indicator?

An example of a precipitation titration is the determination of chloride ions in a sample using silver nitrate (AgNO3) as the titrant and potassium chromate (K2CrO4) as the indicator.

In this titration, a known volume of silver nitrate solution is slowly added to a solution containing the chloride ions to be determined. When all of the chloride ions have reacted with the silver nitrate to form an insoluble compound (silver chloride), a visible precipitate will form.

Potassium chromate is used as the indicator in this titration because it changes color in the presence of excess silver ions. When the reaction is complete and excess silver nitrate is present, the yellow potassium chromate solution will change to orange. This color change indicates the end point of the titration and indicates that all of the chloride ions have reacted with the silver nitrate.

By accurately measuring the volume of silver nitrate solution required to reach the end point, the chemist can calculate the molarity of chloride ions in the original sample.

does hardness and alkalinity occur if we add salt in water?

Adding salt to water can affect the hardness and alkalinity of the water.

Hardness in water is caused by the presence of dissolved minerals, such as calcium and magnesium, which can react with soap to form soap scum and make water appear hard. Adding salt to water can increase its hardness by dissolving more minerals into the water. The type and amount of salt added will determine the extent to which the water hardness increases.

Alkalinity is a measure of a water's ability to neutralize acids and is caused by the presence of dissolved carbonates, bicarbonates, and hydroxides. Adding salt to water can affect its alkalinity depending on the type of salt used. For example, adding sodium chloride (table salt) to water will not significantly change its alkalinity, but adding sodium carbonate (soda ash) to water will increase its alkalinity.

In general, the effect of adding salt to water on its hardness and alkalinity will depend on the type and amount of salt added, as well as the initial hardness and alkalinity of the water.